



3. Preparation of the implant site.

Compressing the bone with Convex compress screw Ø2/3.3mm.



4. Insertion of the SEVEN implant MF7-13375

The SEVEN implants improves compression during insertion, providing a better initial stability.

Sinus lifting

Ø3 - Ø3.5 - Ø4



Instrument Maintenance

Disinfection

- Immerse instruments immediately after use.
- Use approved agents only.
- Observe manufacturer's recommendations regarding concentration / time / material compatibility.

Cleaning

- Remove all residue.
- Use Ultrasound.
- Use anticorrosive cleaning agent.
- Thoroughly rinse cleaning and disinfecting agents under running water.
- Use distilled water to prevent water spots.

Drying

- Dry only with:
 - Compressed air.
 - Hot air.
 - Absorbent paper tissue.

Examination

- Perform visual inspection.
- Dispose of damaged instruments.

Check for:

- Breakouts in blades.
- Bent instruments.
- Corrosion.

Sterilization

- All dental instruments can be sterilized.
- The instruments must be sterilize before use.
- The device must be sterilized before use by autoclave, at a temperature of 134°C (273°F) at a pressure of ≈315 Kpa during 6 minutes.
- Do not exceed 134°C during sterilization.

Storage

- Store in dry, dust-proof area.
- Keep instruments separated from chemicals.

⚠ Attention, see instructions for use

LOT Batch Code

MW Date of Manufacture

REF Catalogue Number

Manufacturer

Non-Sterile

MP-UI008 Rev.6



MIS Implants Technologies Ltd.
P.O.Box 7 Bar Lev Industrial Park 20156, ISRAEL
Website: www.mis-implants.com

MIS®

CE
0483

EC REF

MIS Implants Technologies GmbH
Paulinenstr 12A, 32427 Minden, Germany
Tel: +49 571-972-7690
Email: service@mis-implants.de

© MIS Corporation. All rights Reserved.

User Manual
Bone Compression Kit



MIS®

*Example using SEVEN implant MF7-13375

The Compression bone kit has been designed for making the cavity used by dental surgeon during the operation of Dental Implant, which is benefit for increasing in bone density for soft bone.

- Atraumatic osteotomy.
- Increase in bone density.
- Increase in a primary stability.
- Gradual control of expansion.

Indication for use:

- For internal sinus floor elevation procedure in case the bone level in the lateral maxilla is between 7-10 mm.
- Using the osteotoms is depending the bone type.
- Control of the osteotom –screw hand reduces a minimum risks of perforation in the upper maxilla.
- Laser marks on the screw indicate a height of 8,10,13 mm.



Instructions for use


1.
By hand

2.
By ratchet

3.
By motor

- MO-CO160 | Convex compress screw 1.6/2.8 mm
- MO-CO200 | Convex compress screw 2/3.3 mm
- MO-CO240 | Compress screw 2.4/3.6 mm
- MO-CO280 | Convex compress screw 2.8/4.2 mm
- MO-SO300 | Concave sinus screw 3 mm
- MO-SO350 | Concave sinus screw 3.5 mm
- MO-SO400 | Concave sinus screw 4 mm
- MT-MMA10 | Short motor mount adapter
- MT-MMA20 | Long motor mount adapter
- MT-MRL10 | Long ratchet or hand wrench adapter
- MT-MRS10 | Short ratchet or hand wrench adapter
- MT-SMD10 | Spade marking drill
- MT-RI030 | Ratchet wrench

Surgical sequence



Flowchart

Compress	ØApical	Length			ØImplant
		8	10	13	
1	1.6	Ø=2.3	Ø=2.5	Ø=2.8	3.3
2	2.0	Ø=2.5	Ø=2.7	Ø=3.3	3.75
3	2.4	Ø=3.1	Ø=3.3	Ø=3.6	4.2
4	2.8	Ø=3.5	Ø=3.8	Ø=4.2	4.7-5


*Dimension (mm)

Step by step procedure.

1.
Preparation of the implant site.
Drilling with pilot drill Ø2mm.



2.
Preparation of the implant site.
Compressing the bone with Convex compress screw Ø1.6/2.8mm.



•Example using SEVEN implant MF7-13375