

Original
MD-CTP10-5
MD-CTP10
 Narrow cementing post

Original
MD-MAC10-6
MD-MAC10
 Standard cementing post

Original
MD-CPK81-1
MD-CPK81
 Anatomic 1mm transgingival abutment, h.8mm, S,P

Original
MD-CPK82-2
MD-CPK82
 Anatomic 2mm transgingival abutment, h.8mm, S,P

Original
MD-CPK83-3
MD-CPK83
 Anatomic 3mm transgingival abutment, h.8mm, S,P

Original
MD-CPK84-4
MD-CPK84
 Anatomic 4mm transgingival abutment, h.8mm, S,P

Original
MD-A0010-8
MD-A0010
 Esthetic abutment internal hex. 1mm

Original
MD-AN251-10
MD-AN251
 25° Anatomic angulated abutment

Original
MD-AN151-9
MD-AN151
 15° Anatomic angulated abutment

Original
MD-P1530-13
MD-P1530
 Esthetic 15° angulated abutment internal hex. 3mm

Original
MD-A2510-12
MD-A2510
 Esthetic 25° angulated abutment internal hex. 1mm

Original
MW-CPK81-15
MW-CPK81
 Anatomic 1mm transgingival abutment, h.8mm, W,P

Original
MW-CPK82-16
MW-CPK82
 Anatomic 2mm transgingival abutment, h.8mm, W,P

Original
MW-CTP10-19
MW-CTP10
 Direct conical titanium post for W,P

Original
MW-MAC10-20
MW-MAC10
 cementing post W,P

Original
MW-AN151-21
MW-AN151
 Wide 15° anatomic angulated abutment

Original
MW-P1510-22
MW-P1510
 Wide esthetic 15° angulated abutment internal hex.1mm

Original
MW-CPK83-17
MW-CPK83
 Anatomic 3mm transgingival abutment h.8mm, W,P

Original
MW-CPK84-18
MW-CPK84
 Anatomic 4mm transgingival abutment h.8mm, W,P


Original
MD-P0030-7
MD-P0030
 Esthetic abutment internal hex. 3mm

Original
MD-A1510-11
MD-A1510
 Esthetic 15° angulated abutment internal hex. 1mm

Original
MD-P2530-14
MD-P2530
 Esthetic 25° angulated abutment internal hex. 3mm

Original
MD-GT190
 Guide pin for technician set

* To order items from the planing kit, please **indicate the colored catalog number** in the wall chart.




1a,1b. Indicate the right angulation of implant and adjustment of the abutment by the planning guide pin.

2. Secure prosthetic planning abutment (replica) to implant with planning guide pin.





Planning kit abutments are replicas only.

Prosthetic Planning Kit MK-0040



Standard

Wide



Key to code used:

Non sterile

LOT Batch Code

REF Catalogue Number

Manufacturer



MC-PPK10C Rev. 1



MIS Implants Technologies Ltd.
P.O.Box 7, Bar Lev Industrial Park 20156, ISRAEL
Website: www.mis-implants.com



EC REP

MIS Implants Technologies GmbH
Paulinenstr. 12 A, 32427 Minden, Germany
Tel: +49 571-972-7690
Email: service@mis-implants.de