










DO NOT EXCEED 134°C (273°F) DURING STERILIZATION

mis

Prosthetic Tools Kit Contents:

- 1  **MT-HHR13**
Long hand screwdriver for 0.05" hex.
- 2  **MM-LSD15**
Locator torque wrench driver short
- 3  **MM-LSD21**
Locator torque wrench driver long
- 4  **MT-MUK02**
Multi-unit key
- 5  **MT-IT100**
Closed tray impression coping key
- 6  **MT-BK235**
Ball attachment anchor key
- 7  **MT-HDL30**
Long hex. drive 0.05"
- 8  X 4
MT-DI001
Implant Direction Indicator
- 9  **MT-UW002**
Wide int. hex. Implant analog for universal hand driver
- 10  **MM-LCT10**
Locator core tool
- 11  **MT-RI040**
Torque wrench
- 12  **MT-MSD20**
Surgical Adapter Driver

MT-HHR07
Short hand screwdriver for 0.05" hex.

MT-LMK13
Long motor key adapter 0.05"

MT-SMK13
Short motor key adapter 0.05"

MT-IT100
Closed tray impression coping key

MT-BK235
Ball attachment anchor key

MT-HDL30
Long hex. drive 0.05"

MT-HDS30
Short hex. drive 0.05"

MT-DI001
Implant Direction Indicator

MT-US002
Int. hex. implant analog for universal hand driver

MM-LCT10
Locator core tool

MT-RI040
Torque wrench

MT-MSD20
Surgical Adapter Driver

Key to code used:



 Non sterile

 Manufacturer

!

The devices must be sterilized before use by autoclave, at a temperature of 134°C (273°F) at a pressure of ≈315 Kpa during 6 minutes.

 MIS Implants Technologies Ltd.
P.O.Box 7, Bar Lev Industrial Park 20156, ISRAEL
Website: www.mis-implants.com

 
MIS Implants Technologies GmbH
Paulinenstr. 12 A, 32427 Minden, Germany
Tel: +49 571 972-7830
Email: service@mis-implants.de

MIS[®]

MIS Implants Technologies Ltd.
www.mis-implants.com

MIS's Quality System complies with international quality standards: ISO 13485:2003 - Quality Management System for Medical Devices, ISO 9001:2008 - Quality Management System and CE Directive for Medical Devices 93/42/EEC. MIS's products are cleared for marketing in the USA and CE approved.

© MIS Corporation. All Rights Reserved.

MIS[®]

Prosthetic Tools Kit

wall chart

