

**Crest
Widener Kit.**
MK-CW28

42.

The MIS Crest Widener is used to widen a narrow alveolar crest which has sufficient height for implantation. By separating the cortical plates new bone is created in between the plates, and the attached gingiva is extracted simultaneously. With this technique, bone augmentation is avoided, the implant is placed in the correct lateral position and sufficient attached gingiva is obtained. In addition, the widener can be used to allow easy distal or mesial orthodontic movements of the molars.



The Crest widener consists of 4 arms, 2 on each side connected with guiding pins and an activating screw. Also provided in the kit: handle key and titanium wire.

Advantages.

