

Drill Stopper Kits.

40.

The Drill Stopper Kits ensure the dental surgeon simple and accurate depth control.

Stoppers

Five sleeves for standard implants in 8, 10, 11.5, 13 and 16mm with different-colored rings for each diameter.

IN Insert

Teflon tube simplifies the insertion of the drill into the stopper sleeve.

Drill

Color-coded drill matches sleeve color.

OUT Insert

Teflon tube simplifies the removal of the drill from the stopper sleeve.



Simple.

Quick, easy assembly.

Fast.

Color coding for immediate identification of sleeve diameter.

Safe.

Safe drilling to planned depth.

Clear.

Clearly arranged matching depth sleeves from 8mm to 16mm.

The kit is available in 5 diameters.

Kit Contents:

The MIS Stopper Kit fits MIS drills, enabling the doctor to drill to the exact desired depth during implantation. It is simple to use and available in different diameters (2mm, 2.8mm, 3.2mm, 3.8mm, 4.5mm).



MK-ST020 | Ø 2mm - Standard



MK-ST038 | Ø 3.8mm - Standard



MK-ST028 | Ø 2.8mm - Standard



MK-ST045 | Ø 4.5mm - Standard



MK-ST032 | Ø 3.2mm - Standard